

# Plinth 2000

## Model 502CT Changing Table 1800mm

### Overview:

Designed specifically for use within schools and homes for changing disabled people. Strong, stable platform with hoist access and fitted with fold-down cot sides. Choice of three sizes available: 1200mm and 1500mm with the 1800mm length shown below. Optional padded cot sides available.



### Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Hoist clearance 130mm
- Fold down cot sides

### Overall Dimensions

Length: 180cm (71") approx  
Width: 63cm (25") approx  
Height: 47cm - 97cm approx  
(18" - 38")approx

### Upholstery Dimensions



### Recommended Options:



Padded Cotside




Battery Back-up



Double Foam



Large Castor Base



Paper Roll Holder

NHS Approved Supplier



UK Made

t: +44 (0)1449 767887  
e: sales@plinth2000.com  
www.plinth2000.com

## Operational Parameters

---



Model 502CTE - Electric version

Model 502CTH - Hydraulic version

## For Your Peace of Mind

---

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

## Standards

---

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

## Delivery Promise

---

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

## Upholstery Vinyl

---

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



White



Dove



Black



Chocolate



Damson



Wine



Red



Terracotta



Rose



Night



Ocean



Sky



Slate



Aqua



Green



Sage



Lime



Yellow



Beige

Plinth 2000 Ltd.  
Wetheringsett Manor  
Wetheringsett  
Stowmarket  
Suffolk, IP14 5PP

t: +44 (0)1449 767887  
f: +44 (0)1449 766122  
e: sales@plinth2000.com  
www.plinth2000.com

**Plinth** | 2000