

Plinth | 2000

EXAMINATION & TREATMENT RANGE



Physiotherapy, Osteopathy, Chiropractic, Rehabilitation & Paediatric

CONTENTS:



02	What Counts	11	40Wide Neurology Couch
03	502 2 Section Couch	12	50 Heavy Duty 3 Section Couch
04	503 3 Section Couch	13	501 Tilt Table
05	502 Changing Table	14	502 Traction Table
06	Elite Portable Treatment Table	15	Accessories
07	Kiko Next Generation Couch	17	Upholstery
08	502R Chiropractic Couch		
09	504 5 Section Couch		
10	503M Manipulation Couch		

WHAT COUNTS:



GUARANTEE

Plinth 2000 products have a lifetime guarantee on the frame. All other components (excluding upholstery) are guaranteed for 3 years.



MAINTENANCE

Maintenance is not required on our products although we can arrange regular servicing if you wish.



DELIVERY

Average delivery is just 7 days



PLACE YOUR ORDER

We would love to discuss your order with you to ensure we can meet all your requirements, if however you are unable to phone then just e-mail or fax and we'll do the rest.

Tel: +44 (0) 1449 767 887
Fax: +44 (0) 1449 766 122
Mail: sales@plinth2000.com
Web: www.plinth2000.com

Post: Plinth 2000 Ltd,
Wetheringsett Manor,
Wetheringsett,
Stowmarket, Suffolk,
IP14 5PP, England

502 2 Section Plinth

Robust Two Section couch ideally suited for continuous use in areas such as outpatient departments, gp surgeries, colleges, private clinics etc.

Features: unique lifetime guarantee, easy operation, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

503 3 Section Plinth

Very popular and versatile three-section couch, firm and stable providing the ideal base for general treatment techniques through to heavy manipulation. Commonly used in areas such as physiotherapy, gp surgeries, alternative therapies, private clinics etc. A wide range of accessories can be fitted to expand the versatility of this couch.


Features: unique lifetime guarantee, easy operation, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.





502 DIMENSIONS

Length: 188 cm (74")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

502 Electric
 502 Hydraulic

Model shown with optional breathing hole.
 Please specify motorised back rest or other optional extras if required



-  225kg SAFE WORKING LOAD
-  ANTI MICROBIAL VINYL
-  LEASING AVAILABLE
-  TWO MAN DELIVERY/INSTALLATION





503 DIMENSIONS

Length: 185 cm (73")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

503 Electric
 503 Hydraulic

Model shown with optional breathing hole.
 Please specify motorised back rest or other optional extras if required



-  225kg SAFE WORKING LOAD
-  ANTI MICROBIAL VINYL
-  LEASING AVAILABLE
-  TWO MAN DELIVERY/INSTALLATION



Call +44 (0) 1449 767 887

502 Changing Table

ELITE Portable Treatment Table

Designed specifically for use within schools and homes for changing disabled people. Strong, stable platform fitted with fold-down cot sides. Optional padded cot sides available. Choice of lengths 1200, 1500 & 1800.

Features: cot sides included, hoist accessibility, unique lifetime guarantee, easy operation, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

Elegant portable couch manufactured to the highest standards, fully CE marked and TUV tested. Weighing in at only 16kg this couch is able to offer a maximum user weight of an amazing 160kg!


Features: Thick comfortable upholstery, strong lightweight aluminium frame, weighs only 16kg, height adjustable from 63cm to 76cm, able to hold 160kg, fast & easy folding/unfolding, breathing hole included.

502 DIMENSIONS

Length: 180 cm (71")
 150 cm (59")
 120 cm (47")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

502CT Electric
 502CT Hydraulic

Please specify length when ordering.
 Optional padded Cot sides available.




- 225kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/ INSTALLATION



ELITE DIMENSIONS

Length: 185 cm (73")
 Width: 63 cm (25")
 Height: 63 cm (25") approx
 to 76 cm (30") approx

Colours Available: Dark Blue only



- 160kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE



Call +44 (0) 1449 767 887

KIKO Next Generation 3 Section Couch

502R Chiropractic Couch

This new model combines contemporary styling and versatile design to offer the very strongest of couches able to withstand very heavy manipulation with no flex or whip. It is far stronger and lighter than any other couch available and has super smooth gas cylinder operation and 'soft touch' upholstery for maximum comfort.

Features: unique lifetime guarantee, easy operation, contemporary styling, exceptionally strong frame design, low running costs, 'soft touch' fire retardant foam, able to lift 225kg from its lowest height, full choice of accessories and colours.

New product designed in conjunction with Chiropractors, this table offers a high quality, economic, robust design guaranteeing many years of trouble free operation and maximum patient comfort.


Features: smooth electric operation, firm comfortable upholstery, adjustable foot support, unique lifetime guarantee, easy operation, low running costs, high density fire retardant foam, hardwearing vinyl, able to lift 225kg from its lowest height, full choice of colours.




KIKO DIMENSIONS

Length: 194 cm (76")
 Width: 68 cm (27")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

KIKO Electric

Please specify motorised back rest or other optional extras if required



-  225kg SAFE WORKING LOAD
-  ANTI MICROBIAL VINYL
-  LEASING AVAILABLE
-  TWO MAN DELIVERY/INSTALLATION




Call +44 (0) 1449 767 887

502R DIMENSIONS

Length: 189 cm (74")
 Width: 71 cm (28")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

502R Electric 502R Hydraulic

Please specify foot switch or other optional extras if required.



-  225kg SAFE WORKING LOAD
-  ANTI MICROBIAL VINYL
-  LEASING AVAILABLE
-  TWO MAN DELIVERY/INSTALLATION



504 5 Section Couch

Providing additional comfort for prone lying patients, this 5 section is based on our popular Model 503 with the addition of a breathing hole and adjustable arm rests. Ideally suited to manipulation and examination and can have a wide range of accessories fitted.

This couch is also available with drainage facility as model 505.

Features: height adjustable arm rests, unique lifetime guarantee, easy operation, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

503M Manipulation Couch

Profiled upholstery and a comfortable drop arm headrest with thicker padding makes this couch ideally suited to manipulation. The upholstery narrows to the shoulder area providing maximum comfort for the prone lying patient whilst still enabling the practitioner to turn the patient in complete security.


Features: height adjustable arm rests, thicker padding, paper roll holder, profiled upholstery, extending foot slider, unique lifetime guarantee, easy operation, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.

504 DIMENSIONS

Length: 185 cm (73")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

504 Electric
 504 Hydraulic

Please specify foot switch or other optional extras if required.



- 225kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/INSTALLATION




503M DIMENSIONS

Length: 162 cm (64")
 Width: 55 cm (22")
 tapering to 41cm (16")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

503M Electric
 503M Hydraulic

Please specify foot switch or other optional extras if required.



- 225kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/INSTALLATION



Call +44 (0) 1449 767 887

40 Wide Neurology Couch

Developed in conjunction with the Bobath Centre in London, our 'double' plinth provides probably the most firm and stable treatment surface available. Being 102cm wide and lowering to just 40cm it is perfect for use in hospital gyms and physiotherapy department. Now safely lifts from its lowest point 320kg.

Features: unique lifetime guarantee, easy electric operation, fully hoist accessible, low running costs, high density fire retardant foam, breathing hole, hardwearing anti-microbial vinyl, able to lift 320kg from its lowest height, full choice of accessories and colours.

50 Heavy Duty Bariatric Couch

Every department should have at least one heavy duty plinth, able to safely lift a very heavy patient to a suitable treatment height when the need arises. This bariatric couch is suitable for a wide range of uses and can be fitted with arms, divided legs and other accessories to make it into a truly multi-purpose couch. Motorised backrest recommended.


Features: unique lifetime guarantee, easy electric operation, fully hoist accessible, low running costs, high density fire retardant foam, breathing hole, hardwearing anti-microbial vinyl, able to lift 320kg from its lowest height, full choice of accessories and colours.

40 DIMENSIONS

Length: 189 cm (74")
 Width: 102 cm (40")
 Height: 41 cm (16") approx
 to 94 cm (37") approx

40 Electric
 40 Hydraulic

Please specify motorised back rest or other optional extras if required.



- 320kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/INSTALLATION




50 DIMENSIONS

Length: 183 cm (72")
 Width: 79 cm (31")
 Height: 48 cm (19") approx
 to 99 cm (39") approx

50 Electric

Please specify motorised back rest or other optional extras if required.



- 320kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/INSTALLATION



Call +44 (0) 1449 767 887

501 Tilt Table

Lowering to just wheelchair height and featuring electric height and tilt, large toe braked castors, adjustable harnesses and crutch handles, this tilt table also has four button hand control permitting simultaneous control and fine adjustment of both height and tilt functions.

Features: unique lifetime guarantee, easy electrical operation of height and tilt functions, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 160kg from its lowest height, full choice of accessories and colours.

502T Traction Table

Friction free rolling top traction table, ideal for both lumbar and cervical traction (sitting or lying) including all harnesses, masts and a spreader bar. Available with or without an electronic traction unit.

Features: unique lifetime guarantee, easy operation, all masts and harnesses, hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 225kg from its lowest height, full choice of accessories and colours.


501 DIMENSIONS

Length: 185 cm (73")
 Width: 58 cm (23")
 Height: 48 cm (19") approx
 to 97 cm (38") approx
 Tilt angle: Horizontal to almost vertical

501 Electric

Optional extras if required:

Patient Table, Tilting Foot Board, Tilting Patient Table, Adjustable Back Rest.




- 160kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/ INSTALLATION

502T DIMENSIONS

Length: 196 cm (77")
 Width: 63 cm (25")
 Height: 46 cm (18") approx
 to 95 cm (37") approx

Model 502T Electric
 Model 502T Hydraulic
 Model 502T package*

*Including Traction machine and flexion stool



- 225kg SAFE WORKING LOAD
- ANTI MICROBIAL VINYL
- LEASING AVAILABLE
- TWO MAN DELIVERY/ INSTALLATION



Call +44 (0) 1449 767 887

CONTENTS:



02	What Counts	11	40Wide Neurology Couch
03	502 2 Section Couch	12	50 Heavy Duty 3 Section Couch
04	503 3 Section Couch	13	501 Tilt Table
05	502 Changing Table	14	502 Traction Table
06	Elite Portable Treatment Table	15	Accessories
07	Kiko Next Generation Couch	17	Upholstery
08	502R Chiropractic Couch		
09	504 5 Section Couch		
10	503M Manipulation Couch		

WHAT COUNTS:



GUARANTEE

Plinth 2000 products have a lifetime guarantee on the frame. All other components (excluding upholstery) are guaranteed for 3 years.



MAINTENANCE

Maintenance is not required on our products although we can arrange regular servicing if you wish.



DELIVERY

Average delivery is just 7 days



PLACE YOUR ORDER

We would love to discuss your order with you to ensure we can meet all your requirements, if however you are unable to phone then just e-mail or fax and we'll do the rest.

Tel: +44 (0) 1449 767 887
Fax: +44 (0) 1449 766 122
Mail: sales@plinth2000.com
Web: www.plinth2000.com

Post: Plinth 2000 Ltd,
Wetheringsett Manor,
Wetheringsett,
Stowmarket, Suffolk,
IP14 5PP, England


UPHOLSTERY »



All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. Most include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products

Note: Unless specified WADI will be supplied.

NEW COLOURS - BRIGHTEN UP OR STAND OUT



Apple Mint

Lilac

Candy

Cool Blue

Mulled Wine

Cocoa

STANDARD HARD WEARING LEATHER GRAIN - STAIN RESISTANT AND FAST DRYING



Snow

Brume

Nubian

Wadi

Fern

Spruce

NON SLIP - GRIPPY, MORE COMFORTABLE



Memphis

Cornflower

Rainforest

Skye

Daviot

Brechin

ADDITIONAL PRODUCT CATALOGUES AVAILABLE :

CLINICAL
RANGE

THERAPY
& TATTOO
RANGE

SURGICAL
RANGE

PODIATRIC
RANGE

Plinth | 2000

www.plinth2000.com