

Viking™ XL mobile lift

Ideal for bariatric care



With a lifting capacity of 300 kg, the Viking XL is ideal for lifting situations involving patients of size. Its flexible application includes lifting from the floor to horizontal lifting. The Viking range also offers special accessories that are designed for gait training during rehabilitation.

Features and benefits

Two lifting speeds accommodate a variety of different lifting situations

The **Universal TwinBar™ 670 slingbar** offers additional sling space for patients of size

Optional quick release hook system means caregivers can easily and rapidly change a sling bar, add a Liko Stretch for horizontal lifting or install a scale

Optional armrests help **support patients during ambulation** as well as making it easier for caregivers to manoeuvre the lift

Ergonomic hand control for lift operation displays both battery usage and overload indicators

User display highlights device usage, including lift count, service required indicator and overloads

An extensive assortment of accessories is available including slings, slingbars and scales

Viking XL mobile lift



Handles



Optional Universal TwinBar 670 slingbar



Hand control



Control box

Technical Specifications

Maximum Load	300 kg
Weight	39.9 kg
Height	1500 to 2110 mm
Protection class	IPX4
Slingbar options	Universal TwinBar 670 optional Crossbars
Base	
height	140 mm
length	1420 mm
width	815 to 1130 mm
expansion	electric
Castors	
front	100 mm double band
rear	125 mm double band
Material	aluminium
Emergency lifting	electric
Emergency lowering	manual and electric
Lifting speed (without load)	17 mm / sec or 23 mm / sec
Lifting interval	1370 mm
Maximum lifting height	1910 mm

Sound level	51 dB(A)
Battery	lead acid gel or lithium ion
Charger	built in
Hand control	wired
Protection class hand control	IPX4
Product number	2040043

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Viking XL is a versatile lift intended mainly for use in health care, intensive care and rehabilitation. Class I
Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN183313-02, 12 May 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom