

SARA STEDY ENHANCED STANDING AID



A VERSATILE STAND AND TRANSFER AID

Sara Stedy®, a stand aid that provides mobility-promoting support and encourages more mobile patients and residents to stand up independently.

With an increased inner width, a new seat shape and a safe working load raised to 182 kg, *Sara Stedy* also meets the standing assistance needs of larger

patients. Enhanced features including a pivoting divided seat and pedal-operated chassis legs ensure improved comfort and support, and make working routines easier for carers.

...with people in mind

www.ArjoHuntleigh.com

KEY FEATURES



1 New seat design

The innovative pivoting seat improves transfer efficiency and patient stability.



2 Easy-to-open chasis

Pedal-operated chasis legs adjust easily for optimum access during transfers.




3 Built-in seat handgrips

Handgrips integrated in the seat allow the carer to easily turn the support aid.



4 Better knee support

Provides extra stability or firm support when a patient cannot fully control leg movements.



5 Total-lock castors

New total-lock castors provide complete stability when the aid is stationary.



6 Lower chasis

Smaller castors and redesigned legs mean a lower chasis for better access.



Product Information	
Weight	29.4 kg (64.82 lbs)
Max. safe working	182 kg (400 lbs)
Total length	905 mm (32 ¾")
Total height	1051 mm (41 ½")
Height of chasis	100 mm (4")

Internal width legs open	729 mm (28 ¾")
Low friction castors, the rear two with brakes	
Lifter – Protection class IPX4	
Approvals	
The product are in compliance with current legislation and applicable product standards	

ArjoHuntleigh is a branch of Arjo Ltd Med. AB.

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2011.

ArjoHuntleigh International Ltd

310-312, Dallow Road Luton, Bedfordshire LU1 1TD, UNITED KINGDOM • Phone: +44 1582 745 800 • Fax: +44 1582 745 866

E-mail: international@arjohuntleigh.com • www.ArjoHuntleigh.com