

Plinth 2000

Model 503TEC Echocardiography Couch

Overview:

This couch has been specifically designed for Echocardiography, with features to enable both excellent positioning for cardiac scanning and comfortable patient seating. The two cut-aways allow easy scanning of the chest wall from both right and left sides of the couch. The height, tilt and back rest are electrically adjusted as standard on the couch which is suitable for both right and left handed sonographers.


Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Electric height & back rest operation
- Trendelenburg tilt
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Hoist clearance 130mm
- Paper roll holder
- Head cushion
- Adjustable handgrip
- Padded side support
- Clinician stool (shaped to fit against couch)
- Integrated foot rest for clinician

Overall Dimensions

Length:	210cm (83") approx
Width:	78cm (31") approx
Height:	46cm - 95cm approx (18" - 37") approx

Upholstery Dimensions




Recommended Options:



Foot Switch



Double Foam



Battery back up



CPR Motor Back


NHS Approved Supplier



UK Made

t: +44 (0)1449 767887
e: sales@plinth2000.com
www.plinth2000.com

Operational Parameters



Trendelenburg Tilt 25°



Model 503TEC - Electric version

For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



Plinth 2000 Ltd.
Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk, IP14 5PP

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com
www.plinth2000.com

Plinth | 2000