

Invacare® | *Bed Accessories*



Complies with
IEC 60601-2-52

Supporting your needs

The *Invacare* bed range features a wide range of accessories which give you the possibility to customise your bed both in terms of design and taste. This allows perfect fittings to personal environments with supporting accessories, tailored to individual client needs.

INVACARE
Yes, you can.®

Bed ends



Victoria bed end

Bed end with removable access panel for **SB™ 755** in beech. Compatible with Line full length side rails.



Susanne bed end

SB 755 beech bed ends, compatible with Line full length side rails.



Anita bed end

Additional wooden bed end design in beech for **SB 755**, compatible with Line full length side rails.

Side rails



Line

Full length aluminium side rail in beech.



Line Extendable

Full length aluminium side rail available in beech. Adjustable length 200-220 cm.




Britt V

Full length wooden side rail available in beech.



Side rail height extender

Height of **Medley Ergo**, Line and **Britt V** full length side rails can be increased by 15 cm.



Medley Ergo wooden side rails

Full length wooden side rails for Invacare **Medley Ergo** range in beech.



Medley Ergo steel side rails

Full length steel side rails for Invacare **Medley Ergo** range.



Verso II side rail

A modern, collapsible steel side rail for use on **Medley Ergo**, **Etude Plus** and **SB 755** beds.

Side rails



Invacare **Scala Basic 2**
Foldable steel side rail 42.7cm height.



Invacare **Scala Basic Plus 2**
Foldable steel side rail 42.7cm height with plastic inserts at the mounting brackets.



Invacare **Scala Medium 2**
Foldable steel side rail 48.3cm height for deeper mattresses.



Invacare **Scala Decubi 2**
Foldable steel side rail 56.2cm height for deep mattresses such as Dynamic systems.

Side rail covers



Cover for Line side rail
Padded covers for Line side rail. Cover remains in place when side rail is lowered. Available in padded, extra padded and extended + 10 cm and + 20 cm.



Net cover for Line side rail
Net covers for Line side rail. Cover remains in place when side rail is lowered. Available as standard and extended + 10 cm and + 20 cm.



Cover for full length side rail
Padded cover for Britt V wooden side rail and **Medley Ergo** metal and wooden side rails. Cover remains in place when side rail is lowered. Available in standard and net.



Cover for Scala metal side rails
Padded cover for three quarter length metal side rails. Cover remains in place when side rail is lowered. Available in standard, extended, and extra padded.



Cover for side rail height extender
Padded cover for side rail height extender placed on Britt V, Line and **Medley Ergo** full length wooden side rails.

Support handles



Support handle

A very stable support handle for getting in and out of bed:

- H25 x W30 cm
- H25 x W80 cm
- H40 x W30 cm
- H40 x W50 cm
- H40 x W95 cm



Swivel support handle


The swivel support handle makes it possible to swing the legs out over the handle before the grip is folded out. Can be used alone or in pairs on the **SB™ 755**.



Support handle for *Octave*®

Makes getting in and out of bed much easier for users.
H40 x W30 cm.

Lifting poles



Lifting pole

Catching handle can be adjusted in both height and depth.



Swivel lifting pole

The lifting pole swings away to the side of the bed. Handle can be adjusted in both height and depth.



Lifting pole

Handle can be adjusted in both height and depth. For use with **Octave** beds.

Other accessories



Distance bumper

Prevents walls from being damaged when the bed is moved. Can be placed on the base frame at the end of the bed or at the side of the bed.



Bumper wheels

The bumper wheel can be mounted horizontally or vertically. Available in three different sizes.



Lowering kit

By using 75 mm wheels, **SB™ 755** can be lowered to 25 cm.



Cover tubes

To be mounted on the bed frame to cover the tubes when the bed **SB 755** mattress platform is extended. Available in 5cm and 10cm.



Limitation crossbar

Prevents the leg support from dipping below horizontal position. Can be placed on the base frame. Available in ivory and aluminium. For **SB 755**.



Extension parts for Etude® Plus

Available in 10 cm or 20 cm extension.



Extension for Octave®

Extension for **Octave** for greater flexibility. 10 cm to 20 cm extension.



Mattress support extension for Medley® Ergo

Extends the length of **Medley Ergo** by additional 15 cm



Wooden adapters – Octave bed

Wooden adapters for mounting **Octave** side rail, comes in beech and cherry.



Side rail lock

Strap for locking the side rail in its lowest position. Can be used with any side rail.



Transportation fitting for Etude Plus


Fittings for easy transport of the **Etude Plus** beds.




Transportation rack

Can be used with **SB 755**. For transporting or storing an assembled bed.


Other accessories



Transport fitting for *Medley Ergo*
For easy storage and transport.



Bed Cradle
Bed cradle can be placed on the bed end to alleviate pressure from the duvet on the legs/feet.



IV Rack
① Free-standing IV rack. Height-adjustable, castors for easy movement.
② Can be placed in the lifting pole tube or on the bed frame.



IV Rack mount fittings
Fittings for mounting the IV rack either in the lifting pole tube or on the bed frame.

Electrical accessories



Hand control
Hand control HB80 and HL80 with green button function. Each function of the HL80 can be individually locked. Available for **SB 755**.



Soft control
Function-based hand control with easy activated buttons. Available for all Invacare beds.



Hand control for *Octave*
Hand control HL85. Each function on the HL85 can be individually locked.



Hand control
Hand control for ***Medley Ergo*** available with or without lockable functions.



Holder for hand control
Flexible holder for hand control, easier access for users. Available for soft control, HB70, HB80, HL80 and HL85 and Dewert hand control. Comes in aluminium and ivory.



T-cable
Can be used when two hand controls are needed.



Battery back-up
To be used in case of power failure*.



Extension cable
2.5 metre extension cable for Linak hand control. Can be used for **SB 755, Etude Plus** and **Octave**.

Electrical accessories



ACP box

Gives healthcare personnel the ability to lock the function of one or more motors. Can be used for **SB™ 755**, **Etude Plus** and **Octave**.



Field power supply

Can be used to adjust the height or any other motor when no power is available. For service personnel only, not to be used with a user in the bed.

Furniture



Bed tray

Stable table/tray, fits directly on the side rails which are full length and wooden (90cm wide).



Pausa

Height-adjustable tilt table with plate that can be adjusted from 72 - 116 cm in height. Max load: 10 kg.



Kauma

Height-adjustable table with plate that can be adjusted from 71 - 113.5 cm in height. Max load: 10 kg.



Sonnet Combi

Bedside table in beech. All horizontal surfaces are in plastic laminate. Practical top shelf, one drawer and a cupboard. The flap can be adjusted from 31 - 93 cm in height.



Back Support

Manual backrest support, with possible head support. Can be adjusted from a 45 - 60 degree angle.

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Bed Accessories - UK - 04/2013

* Please note: Battery back up is for use in emergencies to bring the bed back to a neutral position.

Part numbers and prices are available In Invacare Ltd Price Lists.
For more information, please contact Customer Services on 01656 776222.



Yes, you can.®